
# MedLabQC

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## 1.4 ASIGNAR LOTES DE CONTROL

**Assignment of control materials**

List of registered control materials	
Grandox	15480

**QC materials for GLUCOSA**

Level 1 → Unknown

Level 2 → Unknown

Level 3 → Unknown

**Validation**

Change only the assignment of materials

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**GLUCOSA**  
Entry of QC points  
Date: 05/02/2007 Time: 02:57:18 p.m.  
Unknown: 78 mg/dl  
Unknown: 156 mg/dl  
Unknown: 210 mg/dl  
Operator: lasg  
Comment:   
[OK] [Clear]  
Analytes: GLUCOSA  
Queue : 0 [Exit]

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[OK] [Cancel]  Westgard multiples

# 3. CONFIGURACIÓN

The screenshot shows the MedLabQC software interface. The menu bar includes 'Parameters', 'Configuration', 'File', 'Clipboard', and 'Special functions ?'. A dropdown menu is open under 'Configuration', showing options: 'Backup/restore options' (labeled 'a'), 'Password' (labeled 'b'), 'Statistical risks' (labeled 'c'), and 'Shortcuts' (labeled 'd'). The main area displays three control charts for 'GLUCOSA'. Each chart has a y-axis with a red dot representing a data point. The first chart has a y-axis from 0 to 120, with a red dot at 80. The second chart has a y-axis from 0 to 300, with a red dot at 150. The third chart has a y-axis from 0 to 400, with a red dot at 210. To the right of each chart are 'Mobile limits' and 'Common CV' values. The right panel shows the 'Entry of QC points' form for 'GLUCOSA', with fields for Date (05/02/2007), Time (03:20:48 p.m.), and input fields for 'Unknown' values, 'Operator', and 'Comment'. The bottom toolbar contains icons for 'Data table', 'Daily QC', 'Control limits', 'Print charts', 'Assign lots', 'Archives', 'Backup', 'Queue: 0', and 'Exit'.

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The screenshot displays the MedLabQC software interface. The 'File' menu is open, showing options: 'Diskette restore', 'Export', and 'Import'. Three callout boxes labeled 'a', 'b', and 'c' point to these options respectively. The main window contains three control charts for 'GLUCOSA' with the following data:

Chart	Mobile limits	m	N	Common CV
Top	78,000 mg/dl	78,000 mg/dl	1	Common CV
Middle	156,000 mg/dl	156,000 mg/dl	1	Common CV
Bottom	210,000 mg/dl	210,000 mg/dl	1	Common CV

The right panel shows the 'Entry of QC points' form for 'GLUCOSA' with the following fields:

- Date: 05/02/2007
- Time: 03:20:48 p.m.
- Unknown: [ ] mg/dl
- Unknown: [ ] mg/dl
- Unknown: [ ] mg/dl
- Operator: [ ]
- Comment: [ ]
- Buttons: OK, Clear

The bottom toolbar includes icons for: Data table, Daily QC, Control limits, Print charts, Assign lots, Archives, Backup, Queue: 0, and Exit.

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# 5. CLIPBOARD

The screenshot displays the MedLabQC software interface. On the left, three control charts are shown for different analytes. Each chart has a y-axis with a red mean line, a yellow upper limit, a green lower limit, and a red lower limit. A red dot represents the current data point. The first chart has a mean of 78,000 mg/dl, the second 156,000 mg/dl, and the third 210,000 mg/dl. A 'Clipboard' menu is open over the first chart, showing 'Copy' and 'Paste' options. On the right, the 'GLUCOSA' data entry panel is visible, including fields for Date (05/02/2007), Time (03:20:48 p.m.), and three input fields for 'Unknown' values in mg/dl. The bottom toolbar contains icons for Data table, Daily QC, Control limits, Print charts, Assign lots, Archives, Backup, Queue: 0, and Exit.

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The screenshot displays the MedLabQC software interface. The main window shows three control charts for the analyte 'GLUCOSA'. The top chart has a y-axis from 0 to 120 mg/dl, the middle from 0 to 300 mg/dl, and the bottom from 0 to 400 mg/dl. Each chart has a red mean line, a green mobile limit line, and yellow and red control limit lines. A red dot represents the current data point. A 'Special functions' menu is overlaid on the charts, with callouts 'a', 'b', 'c', and 'd' pointing to specific functions: 'a' points to 'Lock limits', 'b' to 'Fix limits', 'c' to 'Erase/archive levels', and 'd' to 'Clear queue'. To the right of the charts, the 'Mobile limits' and 'Common CV' are displayed for each chart. The bottom right panel shows the 'Entry of QC points' form with fields for Date (05/02/2007) and Time (03:20:48 p.m.), and buttons for 'OK' and 'Clear'. The bottom status bar contains icons for Data table, Daily QC, Control limits, Print charts, Assign lots, Archives, Backup, Queue (0), and Exit.

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**Locking control limits**

From **05/02/2007** to **05/02/2007**

Partial statistics for the interval

#	m	CV	N
C1	76,200	5,4	5
C2	160,750	2,8	4
C3	235,000	7,5	4

**Mobile limits**  
m = 76,200 mg/dl  
CV = 5,4 %  
N = 5  
Common CV

**Mobile limits**  
m = 235,000 mg/dl  
CV = 7,5 %  
N = 4  
Common CV

**GLUCOSA**

**Locking control limits**

From [ ] to [ ]

Partial statistics for the interval

#	m	CV	N
C1			0
C2			0
C3			0

Unknown  
**Mobile limits**

Unknown  
**Mobile limits**

Unknown  
**Mobile limits**

Data table | Daily QC | Control limits | Print charts | Assign lots | Archives | Backup | Queue : 0 | Exit

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# 6.2 FIJAR LÍMITES

The screenshot shows the MedLabQC software interface. A central dialog box titled "Fixing control limits" is open, allowing the user to set control limits for "GLUCOSA". The dialog box contains the following fields and options:

- Unknown:
- Target value:  mg/dl
- Deviations allowed before:
- warning:  action:
- Buttons: Fixed, Mobile


The background shows three control charts for different levels of control. The top chart has a mean of 76,200 mg/dl and CV of 5,4%. The middle chart has a mean of 160 mg/dl and CV of 2,8%. The bottom chart has a mean of 235,000 mg/dl and CV of 7,5%. The x-axis for all charts is labeled "5 5 5 5 5" and "Feb 07".

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OK Cancel Print Edit line Delete all Delete selection Undo

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### History of versions

---

MedLabQC is a major update of qualCont 2.11

Version 3.00 :

- Choice of 1 to 3 control levels for each analyte.
- Each QC point is dated, timed and stored with the initials of the operator.
- Westgard's multirules.
- Analytes sorted in any order.
- When changing the batch of a control material, simultaneous restart of all the relevant analytes.
- Application exchanges through the clipboard or delimited text files.
- Backward compatibility. The data file of the previous version can be imported in MedLabQC.
- Control mode with pre-determined limits ("Fixed limits").
- Possibility to un-validate, correct, and re-validate control points.
- Archiving old batches results in separate files.
- User manual in HTML format.

Version 3.11 :

- Drag and drop for sorting analytes and distributing control materials.
- Queue for data entry (delayed validation).
- Limits definable in probability.

Version 3.21 :

- Common CV method.
- Comma accepted as decimal separator.
- Bug fixed in the names of control materials.

Version 3.22 :

- Bug fixed when importing or pasting with a short date format 'MMM' (month with letter abbreviation).

Version 3.23 :

- Correction of a light error in scales ("common CV" method).

Version 3.24 :

- Check box "Common CV" in the dialog "analytes".
- Paragraph 6-4 in the user manual.

## ReadMe.txt - Notepad

File Edit Format View Help

### Rights and licensing

---

MedLabQC is a freeware, available for distribution by any party. This software is unlimited and free to use without licensing. If you enjoy it, please forward a copy to your friends or colleagues.

You may not alter this software in any way, including changing or removing any messages, windows, or files in the release package. You may not decompile, disassemble or reverse-engineer this software.

### Platform

---

windows 95/98 or NT4. The first versions of win95 and NT4 have been released with a library file ComCtl32.dll 4.70 which needs to be updated to 4.72 for working with MedLabQC.

Screen resolution :

- minimum 800x600 pixels
- screen font (96 dpi)

### Distributed files

---

Downloaded file : MedLabQCFr.exe or MedLabQCUK.exe  
(English and French versions).

Self-extraction of these files creates a sub-directory :

C:\MedLabQC\

in which are exploded the following files :

- MedLabQC.exe : Programme.
- ReadMe.txt : This explanation file.
- \*.htm and \*.gif : 4 files for the user manual.

### Installation

---

You must launch C:\MedLabQC\MedLabQC.exe to start the programme. A setup option is put forward on the first opening. It creates a shortcut on the desktop or in the start menu.

### Uninstallation

---

The only modification made to your computer is the creation of the directory where the programme is installed (default C:\MedLabQC\). No DLL file is added nor any entry in the windows registry. To uninstall MedLabQC simply erase the subdirectory where it is installed and the shortcuts that you could have added on the desktop or in the start menu.

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Copyright (C) 1998-1999-2000 Philippe Marquis